



PRODUCTS IN THE MEDICAL FIELD FOCUSED ON ERGONOMICS

Just think, if the things you use every day could be adapted to your needs. Just think, if they were designed to match your height and way you work. Just think, you would be able to get everything done faster and more easily. You would not have to work in uncomfortable and potentially harmful positions.

This is how we think at Rini. And this is exactly how we design and build our ergonomic products.

**Product catalog 2020
issue 3**



Table Of Contents	Side
1. Operating tables - RiEye Mk2S	3
2. Surgeon chairs - Carl	5
3. Patient chairs & Device tables	7
4. Other MedTec-products	9
5. General Terms and Conditions	15

Welcome to Rini!

We are a Swedish company with innovative research and development and have since mid 1990 been a quality supplier of ergonomic products and solutions to medical professionals. We have our own sales and support organization in Sweden and work mainly with appointed distributors on our export markets.

Our business strategy is simple!

Rini offers ergonomic products and solutions uniquely tailored to different treatment areas. Giving doctors and associated personnel a maximum comfortable and productive work environment.

The close contact we have with many highly qualified users is exciting and a positive challenge to continue to develop our products further - always with customer satisfaction in focus! **Read more at www.rinicompany.com**

Rini is a quality certified supplier according to ISO 9001 and ISO 13485.



Contact us by phone +46 8 594 77170 or E-mail info@rini.se and we will tell you more!

1. Operating tables

Operating table RiEye Mk2S

Product code **Description**

- 100-00105-05 RiEye Mk2S Operating table – R5
- Verical leg support easier for patient to sit in the table straight from the front
 - Smooth and easy for the patient to sit down from the front of the table
 - Leg cushion 400mm long



- 100-00105-06 RiEye Mk2S Operating table – R6
- Headrest with an unique ergonomic design
 - Possible to set the backrest 90° inrelation to the seat
 - Leg cushion 500mm long



- 100-00105-07 RiEye Mk2S Operating table – R7
- All tables has the unique adaptive chassis
 - Can be moved sideways
 - Extra long leg cushion 680mm



Colors selection for standard cushion kits

Product code	Description	R5	R6	R7
140-00202-xx	PERI cushion kit blue	-00	-01	-02
140-00203-xx	PERI cushion kit green	-00	-01	-02
140-00204-xx	PERI cushion kit black	-00	-01	-02
140-00205-xx	PERI cushion kit grey	-00	-01	-02
140-00308-xx	PERI cushion kit yellow	-00	-01	-02
140-00310-xx	PERI cushion kit red	-00	-01	-02
140-00285-xx	PERI extra wide cushion kit black	-00	-01	-02
140-00210-00	PERI special color other than standard			

140-00210-01 **New** ILLUSION cushion kits with upholstery in Polyurethane, PVC and Phthalate free

Accessories for operating table RiEye Mk2S

Accessories marked with * must be factory mounted

Product code	Description
140-00229-01	Electrical headrest *
140-00270-00	Headrest X-Tray, replaces EasyFix (standard) headrest
140-00270-01	Pillow for X-Tray black (300x265mm WxL)
140-00270-02	Pillow for X-Tray grey
140-00270-03	New Headrest S-Tray, replaces EasyFix (standard) headrest
140-00270-04	New Pillow for S-tray (210x180mm WxL)
140-00252-02	New Headrest S-Tray Nippon, replaces EasyFix (standard) headrest
140-00252-05	New Headrest X-Tray FaceDown , replaces EasyFix (standard) headrest
140-00270-06	New Pillow for X-Tray FaceDown
140-00235-02	New Accessory rail for S-Tray and X-Tray
000-03420-16	New Foundation for headrests to make headrest dockable
140-00230-01	Foot control with two joysticks
140-00261-01	Auto access to memory location
140-00266-01	Electrical brake and direction control*
140-00163-00	Anaesthesia armrest, clamp and arm strap included
140-00307-00	Width extension back rest right 410mm
140-00307-01	Width extension back rest left 410mm
140-00232-00	Length extension back rest 120mm
140-00242-01	Wrist support EasyFix
140-00242-02	Wrist support S-Tray and X-Tray
140-00259-00	Driving handle backrest foldable R/L
140-00256-00	Driving handle leg foldable R/L R6 and R7
140-00256-01	New Driving handle leg R/L R5
140-00235-00	Attachment rail for accessories on the seat R or L 250mm
140-00235-01	Attachment rail for accessories on the seat R or L 250mm if Side rail is mounted
140-00236-00	Drape support
140-00236-01	Drape support, formable with air tube connection
140-00237-02	Extra battery
140-00238-01	Holder for hand control unit flexible, bracket included
140-00238-02	Holder for hand control unit fixed, bracket included
140-00238-03	Holder for hand control unit fixed for IV pole, bracket included
140-00239-01	IV pole, bracket included
140-00240-00	Side rail collapsible right side
140-00240-01	Side rail collapsible left side
140-00320-00	New Armrest for Side rail collapsible
140-00240-02	New Side rail right or left side 600mm, 140-00235-00 included
140-00251-00	Leg cushion cover R6 R7
140-00251-01	Leg cushion cover R5
140-00263-00	Belt waist
140-00263-01	Belt arm 700mm
140-00263-02	Belt leg 2100mm
140-00265-00	Leg support circular 150mm
140-00265-01	Leg support half circular 200mm

2. Surgeons chairs

Surgeon chairs Carl MK2, a new generation of electrically controlled chairs

Product code **Description**

- 120-00133-03 **New** Carl Mk2 R5 - Easy arm support
- Uniquely designed frame with three wheels for high availability
 - Electric controlled central brake
 - Seat height from 490 to 800mm
 - Seat 420x460mm back support 300x350mm (LxW)
 - Battery-operated, charger is connected to chair
 - Back and seat can be adjusted in height and angle
 - The chair is equipped with Easy arm support



Product code **Description**

- 120-00133-00 Carl Mk2 R6 - Rilis arm support
- Uniquely designed frame with three wheels for high availability
 - Electric controlled central brake
 - Seat height from 490 to 800mm
 - Seat 420x460mm back support 300x350mm
 - Battery-operated, charger is connected to chair
 - Back and seat can be adjusted in height and angle
 - The chair is equipped with Rilis safety arm support



Product code **Description**

- 120-00133-01 **New** Carl Mk2 R7 - One hand grip support
- Uniquely designed frame with three wheels for high availability
 - Electric controlled central brake
 - Seat height from 490 to 800mm
 - Seat 420x460mm back support 300x350mm
 - Battery-operated, charger is connected to chair
 - Back and seat can be adjusted in height and angle
 - The chair is equipped with one hand grip support



Product code **Description**

- 120-00133-02 Carl Mk2 R8 - FlexiDoc
- Uniquely designed frame with three wheels for high availability
 - Electric controlled central brake
 - Seat height from 490 to 800mm
 - Seat 420x460mm back support 300x350mm
 - Battery-operated, charger is connected to chair
 - Back and seat can be adjusted in height and angle
 - With RoboDoc armrests for the most demanding procedures



Upholstery and colour (state when ordered)

Product code	Description	Blue	Black
140-00118-xx	Cushion kit vinyl	-00	-02
140-00247-xx	Cushion kit art leather	-00	-02
140-00306-xx	Cushion kit soft padding	Not offered	-02
140-00321-xx	New Saddle seat Amazone	Not offered	-02
140-00322-xx	New Saddle seat Jumper	Not offered	-02

Vinyl has an airy, ventilated slick texture and is not so hot, good for long sessions
Artificial leather has a denser structure and are easier to clean

Options for One hand grip armrest R7

Product code **Description**

- 140-00323-00 **New** Arm cushion soft and convex shape
140-00324-00 **New** Extended height adjustment 350mm (225mm)

Option for Carl Mk2

Product code **Description**

- 140-00327-00 **New** Pedal platform 600x680mm (WxL) can be retrofitted

Surgeon chairs Carl, electrically and manually controlled chairs

- | | |
|---------------------|--|
| Product code | Description |
| 120-00102-00 | Carl Foot R6 - foot controlled electrical |
| 120-00102-01 | Carl Heel R6 - heel controlled electrical |
| 120-00102-02 | New Carl Foot R7 - foot controlled electrical |
| 120-00102-03 | New Carl Heel R7 - heel controlled electrical |
- Low friction, large 100mm castors
 - Seat height from 500 to 800mm
 - Seat 420x460mm back support 300x350mm (LxW)
 - Back and seat can be adjusted in height and angle
 - Battery-operated, charger is connected to chair
 - Foot controlled brake, brakes at both rear wheels



Foot R6



Heel R7

- | | |
|--------------|--|
| 120-00100-00 | Carl Spring R6 - foot controlled with H-chassis |
| 120-00100-01 | New Carl Spring R7 - foot controlled with H-chassis |
- Low friction, large 100mm castors
 - Foot controlled height adjustment with ring
 - Seat height from 540 to 680mm
 - Seat 420x460mm back support 300x350mm
 - Seat and back can be adjusted in height and angle
 - Foot controlled brake, brakes at both rear wheels



R6



R7

- | | |
|--------------|--|
| 120-00101-00 | Carl Swing R6 - foot controlled and rotatable |
| 120-00101-01 | New Carl Swing R7 - foot controlled and rotatable |
- Foot controlled height adjustment with ring
 - Seat height from 500 to 700mm
 - Seat 420x460mm back support 300x350mm
 - An easily moved and rotatable chair
 - Back and seat can be adjusted in height and angle



R6



R7

Upholstery and colour (state when ordered)

Product code	Description	Blue	Black
140-00118-xx	Cushion kit vinyl	-00	-02
140-00247-xx	Cushion kit art leather	-00	-02

Vinyl has an airy, ventilated slick texture and is not so hot, good for long sessions
Artificial leather has a denser structure and are easier to clean

3. Patient chairs & Device tables

Patient chairs Malin - electric controlled

Product code Description

- 120-00103-00 Malin Rock - hand controlled
- Low friction large 100mm castors with double ball bearing wheels
 - Lifting capacity from 500 to 800mm
 - Seat 450x470mm back support 370x400mm (LxW)
 - Battery powered
 - Stable manage patient weights up to 120Kg



- 120-00103-01 Malin Swing - hand controlled
- Seat can be rotated $\pm 90^\circ$
 - Low friction large 100mm castors with double ball bearing wheels
 - Lifting capacity from 500 to 800mm
 - Seat 450x470mm back support 370x400mm
 - Battery powered
 - Stable manage patient weights up to 120Kg



- 120-00103-02 Malin Slim - hand controlled
- Wheelbase for narrow floor space under the table
 - Low friction large 100mm castors with double ball bearing wheels
 - Lifting capacity from 500 to 800mm
 - Seat 450x470mm back support 370x400mm
 - Battery powered
 - Stable manage patient weights up to 120Kg



Options:

Product code Description

- 140-00128-00 Extended lifting capacity from 540mm to 940mm
 140-00126-00 Wire less hand control infrared 1ch
 140-00126-01 Wire less hand control infrared 2ch (if electrical brake)
 140-00250-00 Electrical brake by the standard hand control
 140-00315-00 Wheels with direction control

Accessories:

Product code Description

- 150-00102-00 Armrest standard left
 150-00102-01 Armrest standard right
 150-00103-00 Armrest collapsible left
 150-00103-01 Armrest collapsible right
 150-00105-00 Armrest for fluorescein photography left 350mm
 150-00105-01 Armrest for fluorescein photography right 350mm

- 140-00129-00 Foot support plate of rubber, not for Malin Slim
 140-00130-00 Foot support T-shaped of metal, not for Malin Slim
 140-00132-00 Driving handle 450mm
 140-00132-02 Driving handle 600mm
 140-00260-00 Belt for patient
 190-00018-00 **New** Head support



Colours of cushion kit:

Product code Description

- 140-00135-00 Cushion kit blue
 140-00135-01 Cushion kit green
 140-00135-02 Cushion kit black

Device table Linnea

Product code Description

- 130-00102-00 Linnea 650 - examination table
- Width between columns 650mm
 - Table top of beech laminate, size 800x350x22mm
 - Electrically adjustable height from 690 to 1040mm
 - Manage weight up to 130Kg
 - Connected to 230 V
 - Equipped with safety brake for the patient

- 130-00103-00 Linnea 850 - examination table
- Width between columns 850mm
 - Table top of beech laminate, size 1000x420x22mm
 - Electrically adjustable height from 690 to 1040mm
 - Manage weight up to 130Kg
 - Connected to 230 V
 - Equipped with safety brake for the patient

- 130-00104-00 Linnea 1050 - examination table
- Width between columns 1050mm
 - Table top of light grey laminate, size 1180x350x22mm
 - Electrically adjustable height from 690 to 1040mm
 - Manage weight up to 130Kg
 - Connected to 230 V
 - Equipped with safety brake for the patient

- 130-00109-00 Linnea Piccolo - examination table,
- Width between columns 720mm
 - Table top of light grey laminate, size 870x365/420x22mm
 - Electrically adjustable height from 580 to 1240mm
 - Manage weight up to 80Kg
 - Connected to 230 V
 - Equipped with low wheel base and two sockets 230V outlets as standard



Examination table Linnea 650 table top variants:

Product code Description

- 140-00140-00 Table top of beech laminate, 800x450x22mm
 140-00141-00 Table top of beech laminate, 1000x420x22mm, asymmetric

Examination table Linnea 850 table top variants:

Product code Description

- 140-00142-00 Table top of beech laminate 1000x350x22mm
 140-00143-00 Table top of beech laminate 1000x500x22mm
 140-00144-00 Table top of beech laminate 1180x350x22mm, asymmetric
 140-00145-00 Table top of grey laminate with birch edge 1000x450x38mm

Examination table Linnea 1050 table top variants:

Product code Description

- 140-00146-00 Table top of light grey laminate, size 1180x420x22mm

Examination table Linnea Piccilo table top variants:


Product code Description

- 140-00150-00 Table top of light grey laminate, size 870x420x22mm
 140-00145-03 **New** Table top of grey laminat with birch edge 1000x420x38mm cut out for H-S chin support

Options:

Product code Description

- 140-00139-00 Low wheel base, table lowered with 75mm
 140-00121-00 Control box mountet to column for more space beneath table top
 140-00248-00 Factory mounted power outlet under table top



Accessories:

Product code Description

- 140-00272-00 Chin support bracket H-S 18mm
 140-00272-01 Chin support bracket H-S 60mm

4. Other MedTech-products

Special chairs for surgeons and assistants

Product code	Description
120-00115-01	ComfortMove - foot controlled
120-00115-00	ComfortMove - hand controlled
	<ul style="list-style-type: none"> • Ideal for long seatings, popular among anaesthetists • Thick and comfortable padding • Seat 430x400mm back support 300x290mm (LxW)

Standard castors

Castors 50mm low friction

Gas spring	Low	Medium	High
Max seat height,	570	700	820
Min seat height,	430	500	580



120-00105-01	RiFlex - foot controlled
120-00105-00	RiFlex - hand controlled
	<ul style="list-style-type: none"> • Adjustable lower back support • Popular Assistant chair for eye or ear surgery • Seat 360x390mm back support 230x370mm

Standard castors

Castors 50mm low friction

Gas spring	Low	Medium	High
Max seat height,	570	690	870
Min seat height,	430	500	580



120-00104-01	Dynamic - foot controlled
120-00104-00	Dynamic - hand controlled
	<ul style="list-style-type: none"> • Chair with unique dynamic lower back support • Popular operator chair within dental care • Seat 400x440mm back support 160x140mm

Standard castors

Castors 50mm low friction, brake unloaded

Gas spring	Low	Medium	High
Max seat height,	600	720	900
Min seat height,	460	530	620



120-00111-01	RiRest - foot controlled
120-00111-00	RiRest - hand controlled
120-00111-03	RiRest - backrest foot controlled
120-00111-02	RiRest - backrest hand controlled
	<ul style="list-style-type: none"> • Soft seat small and movable • Popular within operating theaters • Seat ø330mm back support 230x370mm

Standard castors

Castors 50mm high friction

Gas spring	Low	Medium	High
Max seat height,	600	720	900
Min seat height,	460	530	620



120-00111-01

120-00111-02

Upholstery

Description	Stamskin	Artleather	Leather	Vinyl
ComfortMove	Option	Std	Option	Option
CorrectSit	Std	Option	Option	Option
RiFlex	-	-	-	Std
Dynamic	Std	Option	Option	Option
RiRest	-	-	-	Std

Colour

Product code	Description	Blue	Green	Black	Other colour
	ComfortMove	Option	Option	Std	Option
	CorrectSit	Option	Option	std	Option
140-00245-xx	RiFlex	-00	-01	-02	-
	Dynamic	Option	Option	Std	Option
140-00246-xx	RiRest	-00	-01	-02	-

Choice of castors

Product code	Description
140-00122-00	Castors 50mm low friction - set
140-00122-01	Castors 50mm high friction - set

Choice of gas spring

Product code	Description
140-00124-00	Gas spring low hand controlled chair
140-00124-01	Gas spring medium hand controlled chair
140-00124-02	Gas spring high hand controlled chair
140-00124-03	Gas spring low foot controlled chair
140-00124-04	Gas spring medium foot controlled chair
140-00124-05	Gas spring high foot controlled chair

Accessories

Product code	Description
140-00183-00	Foot support foot controlled chairs
140-00182-00	Foot support hand controlled chairs

Saddle chairs for surgeons and assistants

Product code Description

120-00106-01	Amazone - saddle chair foot controlled
120-00106-00	Amazone - saddle chair hand controlled
120-00106-02	Amazone - back support hand controlled
120-00106-03	Amazone - back support foot controlled
120-00106-04	AmazoneBalance - hand controlled
120-00106-05	AmazoneBalance - foot controlled

- For you who find other saddle chairs too wide
- Thick and comfortable padding and tiltable seat
- Seat 340x350mm back support 210x330mm (LxW)

Standard castors

Castors 65mm break unloaded

Gas spring	Low	Medium	High
Max seat height,	570	690	820
Min seat height,	430	500	580



120-00106-03



120-00106-05

120-00110-01	Jumper - saddle chair foot controlled
120-00110-00	Jumper - saddle chair hand controlled
120-00106-02	Jumper - back support hand controlled
120-00106-03	Jumper - back support foot controlled
120-00106-04	JumperBalance - hand controlled
120-00106-05	JumperBalance - foot controlled

- A chair for longer sittings
- Excellent for you who wants a lower seat position
- Seat 360x430mm back support 210x330mm

Standard castors

Castors 65mm break unloaded

Gas spring	Low	Medium	High
Max seat height,	570	690	820
Min seat height,	430	500	580



120-00123-01	Perfect lite - saddle chair foot controlled
120-00123-00	Perfect lite - saddle chair hand controlled

- Easy to maintain clean due to a seamless upholstery
- A small chair, excellent for them who want a high seat position
- Seat 300x390mm

Standard castors

Castors 50mm low friction

Gas spring	Low	Medium	High
Max seat height,	570	690	820
Min seat height,	430	500	580



Upholstery

Description	Stamskin	Artleather	Leather
Amazone	Option	Std	Option
Jumber	Option	Std	Option
Perfect classic	Option	Option	Std
Perfect lite	Option	Option	Spec

Colour

Name	Black	Other colour
Amazone	Std	Option
Jumper	Std	Option
Perfect classic	Std	Option
Perfect lite	Std	Option

Choice of castors

Product code Description

140-00122-00	Castors 50mm low friction - set
140-00122-01	Castors 50mm high friction - set

Choice of gas spring

Product code Description

140-00124-00	Gas spring low hand controlled chair
140-00124-01	Gas spring medium hand controlled chair
140-00124-02	Gas spring high hand controlled chair
140-00124-03	Gas spring low foot controlled chair
140-00124-04	Gas spring medium foot controlled chair
140-00124-05	Gas spring high foot controlled chair

Accessories

Product code Description

140-00183-00	Foot support foot controlled chairs
140-00182-00	Foot support hand controlled chairs

Instrument table Anna - for operating theatres

Product code **Description**

- 130-00105-00 Anna Tripod - instrument table with three wheels
- Height adjustment from 700 to 1100mm
 - Size of table top 700x500mm
 - Individually braking castors, low friction
 - Battery powered, charger connect to table
 - Sound absorbing table top



- 130-00106-00 Anna Cross - instrument table with four wheels
- Height adjustment from 850 to 1350mm
 - Size of table top 1140x600 mm
 - Very stable, can manage up to 75Kg
 - Battery powered, charger connect to table
 - Sound absorbing table top



- 130-00106-01 Anna Cross - instrument table with fall edges and four wheels
- Fall edges so the Medical-guidewire for Endoscopic procedures remains on the table
 - Fall edges designed to prevent bacterial collections
 - Size of table top 1000x600 mm
 - Very stable, can manage up to 75Kg
 - Battery powered, charger connect to table



Instrument table variants:

Product code **Description**

- 140-00149-00 Height adjustment 850-1350mm table 130-00105-00
- 140-00151-00 Table top size 600x400mm table 130-00105-00
- 140-00152-00 Table top size 1300x500mm table 130-00106-00

Armrests for operating theatre and examination

Product code	Description
150-00100-00	Rilis armrest R6 <ul style="list-style-type: none"> • Unique safety armrest • Adjusted with great precision • Brackets fitting all of Rinis operator and assistants chairs



150-00113-00	New One grip armrest R7 <ul style="list-style-type: none"> • Unique adjustment can be adjusted draped • Adjusted with great precision • Brackets fitting all of Rinis operator and assistants chairs
--------------	--



Armrest brackets

Product code	Description
000-03067-00	Bracket right/left elliptic
000-01130-00	Bracket right/left rectangular 30x15mm
140-00303-00	Bracket right/left ComfortMove special chair
140-00303-01	Ny Bracket R/L elliptic with rail
000-06302-00	Ny Clove right side
000-06302-00	Ny Clove left side

Undercarriges

Product code	Description
000-03577-01	Under carriage saddle chairs elliptic <ul style="list-style-type: none"> • Jumper, Classic, MultiAdjuster, Swing, Perfect classic, Perfect lite
000-03864-01	Under carriage Bambach elliptic
140-00200-00	Under carriage RiFlex Bambach
000-01325-01	Under carriage saddle chairs 30x15mm <ul style="list-style-type: none"> • Amazone, Jumper, Perfect classic, Perfect lite
000-03940-01	Under carriage Dynamic elliptic
000-01674-01	Under carriage CorrectSit 30x15mm



Armrest for operation and health care

Arm rests

Product code	Description
000-03379-00	Armrest plate in moulded seamless design
000-01327-00	Armrest plate arched
000-03395-00	Armrest plate 270x130mm
000-00431-00	Armrest plate 345x170mm



Arm

Product code	Description
140-00127-00	Arm 75° 400mm, ball joint for curved armrest
140-00180-00	Arm 105° 400mm long



Bracket chair

Product code	Description
000-03374-00	Bracket right elliptic
000-03375-00	Bracket left elliptic
000-03372-00	Bracket right rectangular 30x15mm
000-03373-00	Bracket left rectangular 30x15mm



Armrest Rilis for patient chairs dental

Product code	Description
150-00101-00	Rilis Armrest - for dental unit <ul style="list-style-type: none"> • Unique safety armrest • The armrest plate is upholstered after request • Placement left or right side • Brackets for many different patient chairs



Brackets for dental units

000-03089-00	Bracket flex for dental unit <ul style="list-style-type: none"> • CanCan, Planmeca, HEKA, XO, Stenstrup and more
000-03099-00	Double bracket for KaVo unit
000-06090-00	Bracket for KaVo unit makes the Rilis easy to mount and remove

Ladder Lars

Product code **Description**

- 160-00102-00 Ladder Lars with handle
- Height 300mm at highest step
 - Very stable
 - Provided with anti-slip protection
- 160-00102-01 Ladder Lars without handle
- Height 300mm at highest step
 - Very stable
 - Provided with anti-slip protection
- 160-00102-02 Ladder Lars 3 steps with handle
- Height 450mm at highest step
 - Very stable
 - Provided with anti-slip protection



MEDINOVÀ
> oftalmologia

via Antiniana, 2i | 80078 Pozzuoli | NA (Italy)
tel. +39 081 770 73 27 / +39 081 770 31 98 | fax +39 081 546 32 99
info@medinovasrl.com / www.medinovasrl.com

5. General Terms and Conditions

1. Orders, Prices etc.

All prices are exclusive of VAT and other taxes or official duties. A quotation shall be binding for RINI for thirty (30) days from the date of the offer, unless otherwise stated in writing.

Any information on weight, measurement, capacity etc. in catalogues, prospects, folders are estimates only. All technical data is provided subject to design modifications.

A change of customs, taxes or other official duties, which results in increased cost to RINI, entitles RINI to increase the offered or agreed price.

2. Payment

Orders are to be pre-paid before shipment if not otherwise agreed.

For customers with an established credit line, payment is due strictly thirty (30) days from date of invoice. Any payments not made when due shall accrue interest from the due date until the date payment is made at a per annum rate equal to 18 percentage and will result in that future orders are to be pre-paid.

3. Transportation and Packing

RINI handle packing and delivery of products according to "INCO TERMS Ex Works" or not otherwise agreed.

If other INCO TERMS are agreed for delivery, transportation, insurance and other handling fees will be charged in accordance with RINI's rates prevailing at the time of delivery. If the customer request express delivery, e.g. express cargo, courier, airfreight, the customer is obliged to compensate RINI for the actual cost of transportation in addition to the normal handling charge.

4. Risk of loss of Product

The risk of the product passes to the customer when RINI has delivered the product according to the agreed INCO TERM. Hence, RINI is only responsible for damage that may occur to the product up to the point in time when the risk passes to the customer.

5. Delivery time

The indicated delivery time starts when all terms in the purchase contract has been fulfilled. In case of pre-paid orders the delivery time starts when RINI has received payment. The product is considered to be delivered when risk of loss and damage passes to the customer in accordance with Section 4.

6. Inspection of Product

When a delivery is received the customer shall inspect the transportation packages and products for any transportation damages. If a transport damage has occurred, before the risk of the product has passed to the customer, the customer shall immediately report this to the shipping company and send a copy to RINI with full details of the damages including digital photo.

7. Warranty

RINI warrants that product sold at the time of the passing of risk as stipulated in Section 4 will be free from faults and shortcomings in construction, material and manufacture for 12 months. If the product is of the type requiring alignment or service, RINI is only responsible for faults remaining after such alignment or service has been carried out.

The customer shall immediately after the discovery of the fault or shortcoming make a claim to RINI in the form of a detailed description of the fault or shortcoming. However, such claim must be presented to RINI no later than 7 days after the day the customer discovered or ought to have discovered the fault or shortcoming.

After the receipt of a claim RINI will, at RINI's option, accept the return of the faulty product and reimburse the customer an amount equal to the purchase price, or repair or replace the faulty product.

RINI's responsibility for faults or shortcomings of delivered products is limited to the above. The customer is not entitled to any other remedies, but the above, in relation to faults or shortcomings.

8. Indirect Damages

The customer is not entitled to, and RINI is not in any event liable to pay, compensation for damages which the delivered product has caused to other property or to persons or any other consequential damages, including but not limited to loss of profit, loss of production or any other indirect damages.

In the event RINI would be held liable for compensation to third party for fault or safety shortcomings under the stipulations of, as it is called, product liability, the customer shall fully indemnify RINI for all damages finally awarded or settlement agreed, including RINI's costs in relation thereto, unless the damages caused by the fault or security shortcoming were due to gross negligence of RINI. RINI will give the customer notice of such claim and make it possible for the customer, if customer so wish, to influence the formulation of RINI's statement of defence or the terms to a possible settlement with the person presenting the claim, provided the customer issues a security, acceptable to RINI, for any costs or damages which can be imposed upon RINI.

9. Retention of Title

Products delivered shall remain the property of RINI until payment in full has been made.

10. Intellectual Property Rights

Customer is obliged to immediately notify RINI of any claim from third party alleging that products delivered under this Agreement infringes intellectual property rights (patent, registered design, trademark or copyright) of a third party. Provided customer present such claim immediately and provided customer grants RINI sole authority to defend or settle the case, RINI undertakes, at RINI's expense, to defend the case, or come to a settlement, with the third party whose intellectual property rights has been infringed by the product delivered under this Agreement.

If it is held that a product delivered under this Agreement has infringed third party intellectual property rights RINI shall, at his own expense, and at RINI's option, obtain for the customer the right to continue using the product, replace or modify it so it becomes non-infringing or, in exchange for the return of the product, reimburse the price of the product. RINI is always entitled to deduct and the customer is obliged to pay an amount equal to the use of the product from the time of delivery to the claim was presented. If RINI has modified or developed a new product upon customer instructions or specifications, customer shall fully reimburse RINI for any and all costs and damages that might occur as a result of the modified or new product infringing third party intellectual property rights.

11. Force Majeure

If either party is prevented from performing under this Agreement due to force majeure the other party may not establish any sanctions for this reason. However, if the hindrance regards the customer's possibility to carry out payment RINI is entitled to withhold delivery of the product until payment has been effected. If the hindrance to effect payment last for more than four (4) weeks RINI is entitled to cancel the purchase.

The term "force majeure" means in this Agreement that the performance under the agreement is prevented, aggravated or delayed due to war, rebellion or riot, mobilization or unexpected call for military service of similar dimensions, requisition, seizure, exchange restrictions, export or import restrictions, general scarcity of merchandise, shortage of transportation, labour conflict, restrictions on fuel and power, fire, faulty or delayed deliveries from subcontractors due to above causes, as well as any other circumstance beyond the control of the party.

12. Governing Law and Disputes

This agreement and any disputes hereunder shall be governed by Swedish law.

Any legal action relating to this Agreement and all legal relationships associated with it shall be decided by District Court of Stockholm as exclusive jurisdiction of first instance. However, RINI is entitled to turn to the district court or the Enforcement Service in the county or city where the customer has his legal domicile as regards claim for payment due.

The fact that a dispute has occurred does not exempt the customer to make payment due or to fulfil any other obligations under the agreement.

13. Miscellaneous

Either party's failure to enforce breach of this Agreement or to exercise its rights due to such breach shall not be deemed a waiver of such rights to enforce any future breach of agreement.

What is stated in this document as regards all matters contained is the full and exhaustive description of RINI's responsibilities and the rights of the customer. As regards such matters dealt with the customer is not entitled to any other remedy or compensation.