



For nearly 200 years, Moria has been known for its expertise in making the ophthalmic microsurgical instruments. All our instruments are designed, manufactured and hand finished in France.

The Moria product ranges are covering the following area of expertise:



**Refractive
surgery
equipments**



**Keratoplasty
devices**



**Retina
instruments**



**All surgeries
Steel reusable
instruments**



**All surgeries
Single Use
instruments**



Single-Use Instruments



» SPECULUMS 6



» FORCEPS AND NEEDLE HOLDERS

Capsulorhexis	10
Tissue	13
Monofilament Tying	15
Needle Holders	17



» HOOKS AND RINGS 20



» SCISSORS 26



» CANNULAS AND MICRO-SPONGES 31




» PROCEDURES

Easyset	34
IVT	36
Toric IOL Positioning	37
Neonatology	38
Refractive Surgery	40
Consultation	42



» KERATOPLASTY


Instruments	48
Punches	50
Trephines	52
Artificial Chambers	53
Microkeratome Accessories	54





» SPECULUMS

6




SPECULUMS

► Neonatal Wire Speculum

6 and 9 mm blades

Box of 10 units




17258X10

► Wire Speculum

15 mm blades

Box of 5 units

Box of 10 units



17537X5
17537X10


► Kratz Wire Speculum

15 mm open blades, 1 mm wire

Box of 10 units




17227X10



SPECULUMS

► Adjustable Speculum


15 mm blades
Can be used in nasal and
temporal sides
Box of 10 units




17508X10

► Adjustable Speculum

Solid, 15 mm blades
Box of 10 units




17231X10





» FORCEPS AND NEEDLE HOLDERS

Capsulorhexis	10
Tissue	13
Monofilament Tying	15
Needle Holders	17



CAPSULORHEXIS FORCEPS

► Capsulorhexis forceps - Straight - 90° pointed tips

1.8 mm incision

Box of 10 units



17502X10


► Capsulorhexis forceps - Curved - 90° pointed tips - Corydon type

1.8 mm incision

Box of 10 units



17503X10




CAPSULORHEXIS FORCEPS

▶ Ultrata Forceps - Straight

1.8 mm incision
45° pointed tips
Ergonomic round handle

Box of 5 units
Box of 10 units



NEW




17511X5
17511X10


▶ Crozafor Forceps - Curved

1.8 mm incision
45° pointed tips
Ergonomic round handle

Box of 5 units
Box of 10 units



NEW



17512X5
17512X10


▶ Capsulorhexis forceps - Curved - Corydon type

1.8 mm incision
90° pointed tips
Ergonomic round handle


Box of 5 units
Box of 10 units



NEW



17513X5
17513X10



CAPSULORHEXIS FORCEPS

► Capsulorhexis forceps - Curved - 90° pointed tips

1.8 mm incision

Box of 10 units




17237X10


► Capsulorhexis forceps - Curved - 45° pointed tips

1.8 mm incision

Box of 10 units




17252X10



TISSUE FORCEPS

► Bonn Forceps - Straight


0.12 mm teeth
3.7 mm platforms
Box of 10 units



17221X10

► Bonn Forceps - Straight


0.12 mm teeth
5 mm platforms
Box of 10 units




17504X10

► Bonn Forceps - Angled

0.12 mm teeth
5 mm platforms
Box of 10 units



17505X10




TISSUE FORCEPS

▶ Notched Forceps

Notched jaws 0,3 mm
5mm platforms

Box of 5 units
Box of 10 units



NEW

17514X5
17514X10

▶ Paufigue Forceps

0.4 mm oblique teeth
5 mm platforms

Box of 5 units
Box of 10 units




NEW

17516X5
17516X10

▶ Corneal Graft Preparation Forceps


1.2 mm oval-shaped tip
Curved, 15 mm long

Box of 5 units
Box of 10 units



NEW

17521X5
17521X10





MONOFILAMENT TYING FORCEPS

► Mc Pherson Tying forceps - Angled

8 mm platforms

Box of 10 units




17233X10

► Tying forceps - Straight

4 mm platforms

Box of 10 units





17222X10

► Tying forceps - Curved

4 mm platforms

Box of 10 units




17225X10

3.50 mm platforms

Box of 10 units

17247X10




MONOFILAMENT TYING FORCEPS

► Mc Pherson Tying Forceps - Angled

5 mm platforms

Box of 10 units



17500X10

► Tying Forceps - Curved

5 mm platforms

Box of 10 units



17501X10

► Troutman Forceps - Straight

5 mm platforms


Box of 5 units

Box of 10 units



NEW

17515X5
17515X10




MONOFILAMENT TYING FORCEPS

► Jeweler Forceps - Straight

11.5 cm length

Box of 10 units

NEW



17517X10




NEEDLE HOLDERS

► Needle Holder - Curved


Smooth and curved jaws
Round and ergonomic handle
12 cm length

Box of 10 units

NEW



17534X10





» HOOKS AND RINGS



20

HOOKS AND RINGS

► RELEX manipulator

Curved spatula
Angled hook
Box of 10 units

NEW




17522X10

► LASIK manipulator

Straight and foam tip to move
the corneal flap
Angled hook
Box of 10 units


NEW



17523X10

► Sinskey IOL Manipulator

Angled, 0.2 mm tip
Box of 10 units



17223X10

HOOKS AND RINGS

► Lester IOL Manipulator

Straight, collar button tip

Box of 10 units
FDA only



NEW

17260X10

► Lester IOL Manipulator

Angled, collar button tip

Box of 10 units



17228X10

► Seibel Lens Manipulator

Angled, "Sharp hook" tip shape.
Hook tip angle: 55°

Box of 10 units
FDA only



NEW

17261X10

► Manipulator

Conical, straight, blunt

Box of 10 units




17230X10

HOOKS AND RINGS

▶ Lens Manipulator phaco chopper

1.5 mm tip, inner cutting edge

Box of 10 units




17229X10

▶ Drysdale Lens Manipulator

Blunt paddle design

Box of 10 units




17242X10

▶ Thornton-Fine Fixation Ring

Blunt teeth

1.3 mm internal diameter

Box of 10 units




17232X10

HOOKS AND RINGS

► Price Hook


Reverse Sinskey
Recommended in **DSAEK** and
DMEK for the Descemetorhexis
Box of 5 units



17302X5

► 90° Spatula

90° angled tip
Recommended in **DSAEK** and
DMEK for Descemet's stripping
Box of 5 units



17303X5

► Mendez Ring


- 12 mm inside diameter
- Indicators at increments of 5°
- Four notches on the internal edge at cardinal positions
- Handle positioned at 0° to hold temporally
- Matte surface for maximum visibility under the microscope



Box of 5 units
Box of 10 units




17250X1
17250X10





» SCISSORS

26



SCISSORS

► Scissors - Sharp

Recommended for cutting operative field or removing 6/0 to 10/0 thread

Box of 5 units




17506X5


► Scissors - Blunt

Recommended for cutting operative field or removing 6/0 to 10/0 thread

Box of 5 units



17507X5



SCISSORS

► Vannas Scissors - Curved

Blunt tip
Blade of 5 mm length
Box of 10 units

NEW




17531X10

► Castroviejo Scissors - Curved

Blunt tip
Blade of 8 mm length
Box of 10 units

NEW




17532X10


► Westcott Scissors - Curved

Pointed tip
Blade of 8 mm length
Box of 10 units

NEW



17533X10





» CANNULAS AND MICRO-SPONGES..... 31



CANNULAS AND MICRO-SPONGES

▶ Retrobulbar Needle - 25 G (0.50 mm)

35 mm length, sterile

Box of 10 units




143BIS

▶ Straight Cannula, for children and adults - 25 G (0.50 mm)

Non traumatic and slightly blunt,
12 mm long, outer diameter 25 G,
inner 31 G, sterile

Box of 30 units



13300

▶ Micro-sponges

Absorbent pads with sleeves,
sterile, blunt tip

Box of 20 bags of 10 units



Distributed by:

Moria



30601

▶ Sponcals

Absorbent pads, non-sterile,
cuneiform

Box of 100 units




Distributed by:

Moria




30761






» PROCEDURES

Easyset	34
IVT	36
Toric IOL Positioning	37
Neonatology	38
Refractive Surgery	40
Consultation	42



EASYSET

SPECULUMS	17227	one®	Kratz Wire Speculum
	17508	one® steel	Adjustable Speculum
CAPSULORHEXIS FORCEPS	17237	one®	Capsulorhexis forceps - Curved - 90° pointed tips
	17252	one®	Capsulorhexis forceps - Curved - 45° pointed tips
	17503	one® steel	Capsulorhexis forceps - Curved - 90° pointed tips - Corydon type
	17502	one® steel	Capsulorhexis forceps - Straight - 90° pointed tips
BONN FORCEPS	17221	one®	Bonn Forceps - Straight
	17504	one® steel	Bonn Forceps - Straight
	17505	one® steel	Bonn Forceps - Angled
TYING FORCEPS	17222	one®	Tying forceps - Straight
	17225	one®	Tying forceps - Curved
	17233	one®	Angled Mc Pherson Tying Forceps
	17247	one®	Tying Forceps - Curved
	17501	one® steel	Tying Forceps - Curved
	17500	one® steel	Mc Pherson Tying Forceps - Angled
IOL MANIPULATORS	17223	one®	Sinskey IOL Manipulator
	17228	one®	Lester IOL Manipulator
	17260	one®	Lester IOL Manipulator (FDA only)
LENS MANIPULATORS	17229	one®	Lens Manipulator phaco chopper
	17230	one®	Manipulator
	17242	one®	Drysdale Lens Manipulator
	17261	one®	Seibel Lens Manipulator (FDA only)
SCISSORS	17506	one® steel	Scissors - Sharp
	17507	one® steel	Scissors - Blunt



EASYSET

Cataract sets delivered in sterile blister packs containing a range of 3 to 7 instruments.

EASYSET packs are available in customized sets of both **one**® and **one steel** instruments, containing up to:

- 3 Forceps
- 2 Manipulators
- 1 Speculum
- 1 Scissors



Easy

- One-finger blister pack opening
- Once the blister pack is opened, the instruments are immediately ready for use

Fast

- Fast preparation of the instruments

Convenient

- Compact and ergonomic tray
- Instruments held in their dedicated locations




IVT

► Wire Speculum

15 mm blades

Box of 5 units

Box of 10 units



17537X5
17537X10

► Kratz Wire Speculum

15 mm open blades, 1 mm wire

Box of 10 units



17227X10

► Scleral Marker - 3.5 mm / 4 mm


Scleral marking at 3.5 mm and 4.0 mm

Box of 5 units

Box of 20 units



17539X5
17539X20



TORIC IOL POSITIONING

► Bubble Level

Preoperative marking of the 0-180° axis
Atraumatic and thin marking surface:
0.4 mm x 3 mm
Box of 10 units



17249X10

► Mendez Ring

12 mm inside diameter
Indicators at increments of 5°
Four notches on the internal edge at
cardinal positions
Handle positioned at 0° to hold
temporarily
Matte surface for maximum visibility
under the microscope



Box of 1 unit
Box of 20 units



17250X1
17250X10

► Bores Marker

Used inside the Mendez Ring to make
two radial marks along the desired axis
for precise alignment
Precise 0.4 mm fine markings
Blade design follows curvature of the
cornea
Matte surface for maximum visibility
under the microscope



Box of 10 units



17251X10




NEONATOLOGY

► ROP Kit: Speculum and loop Depressor

Includes a speculum (ref. 17258) and a loop depressor, in sterile pack.

Box of 10 units.



17257X10

► ROP Kit: Speculum and double-ended Depressor

Includes a neonatal speculum and a double-ended depressor (ref. 17253) in sterile pack.

Box of 20 units



17210



NEONATOLOGY

▶ ROP Kit: Speculum and loop Depressor

Includes a neonatal speculum and a loop depressor in sterile pack.

Box of 1 unit
Box of 20 units



17255X1
17255X20

▶ Double-ended Depressor

Box of 20 units




17253X20

▶ Neonatal Wire Speculum

6 and 9 mm blades

Box of 10 units



17258X10



REFRACTIVE SURGERY

► Corneal Marker

Marker for patient's cornea
Available in 8mm, 8,5mm and 9mm diameter
Box of 10 units

NEW



8mm

17518X10

8,5mm

17519X10




9mm

17520X10

► RELEX Manipulator

Curved spatula
Angled hook
Box of 10 units

NEW





17522X10



REFRACTIVE SURGERY

► LASIK Manipulator

Straight and foam tip to move
the corneal flap
Angled hook
Box of 10 units



NEW



17523X10

► Epithelial Spatula

Spatula for the corneal epithelium
Box of 10 units

NEW




17541X10

CONSULTATION

► Tying Forceps - Straight

4 mm platforms

Box of 10 units



17222X10

► Mc Pherson Tying Forceps - Angled

8 mm platforms

Box of 10 units





17233X10

► Tying Forceps - Curved

4 mm platforms

Recommended for removing
10-0 thread or eyelashes

Box of 10 units




17225X10

CONSULTATION

► Mc Pherson Tying Forceps - Angled

5 mm platforms

Box of 10 units



17500X10

► Tying Forceps - Curved

5 mm platforms

Recommended for removing
10-0 thread or eyelashes

Box of 10 units



17501X10

► Troutman Forceps - Straight

5 mm platforms

Box of 5 units

Box of 10 units



NEW

17515X5
17515X10


CONSULTATION

► Bonn Forceps - Straight

0.12 mm teeth, 3.7 mm platforms

Recommended for removing 6/0 or 7/0 thread, grasping compresses or filling cavities with compresses

Box of 10 units




17221X10

► Bonn Forceps - Straight

0.12 mm teeth, 5 mm platforms

Recommended for removing 6/0 or 7/0 thread, grasping compresses or filling cavities with compresses

Box of 10 units




17504X10

► Bonn Forceps - Angled

0.12 mm teeth, 5 mm platforms

Box of 10 units



17505X10

► Scissors - Sharp

Recommended for cutting operative field or removing 6/0 to 10/0 thread

Box of 5 units



17506X5

► Scissors - Blunt

Recommended for cutting operative field or removing 6/0 to 10/0 thread

Box of 5 units




17507X5

► Straight Cannula, for children and adults - 25 G (0.50 mm)

Non traumatic and slightly blunt, 12 mm long, outer diameter 25 G, inner 31 G, sterile

Box of 30 units



13300





» KERATOPLASTY

Instruments	48
Punches	50
Trephines	52
Artificial Chambers	53
Microkeratome Accessories	54



KERATOPLASTY Instruments

► Price Hook

Reverse Sinskey, recommended in **DSAEK** and **DMEK** for the Descemetorhexis

Box of 5 unitss



17302X5

► 90° Spatula

90° angled tip, recommended in **DSAEK** and **DMEK** for Descemet's stripping

Box of 5 units



17303X5

► Busin Spatula

Recommended in **DSAEK**

For insertion of donor graft; 5.5-mm incision

Box of 5 units



17300X5

► 23 G Busin Forceps

Recommended in **DSAEK** for grasping and positioning the donor graft

Box of 1 or 5 units



17301X1
17301X5




KERATOPLASTY Instruments

► Corneal Graft Preparation Forceps

1.2 mm oval-shaped tip
Curved, 15 mm long

Box of 5 units
Box of 10 units



NEW

**17521X5
17521X10**

► Corneal Marker

Marker for patient's cornea
Available in 8mm, 8.5mm and 9mm diameter

Box of 10 units



NEW

8mm

17518X10

8.5mm

17519X10

9mm

17520X10

► Epithelial Spatula

Spatula for the corneal epithelium

Box of 10 units



NEW

17541X10



KERATOPLASTY Punches


► Busin Punch

With aspiration
Straight-walled cut (symmetrical double-bevel blade)
throughout the thickness of the graft


Ideal for DSAEK, DALK, PK

Available in:

- 6.00 mm
- 6.50 – 9.00 mm (increment of 0.25 mm)
- 9.50, 10.00 mm



Blade profile



Shape of the cut


17200DXXX

► Muraine Punch

With aspiration
Straight-walled cut (symmetrical double-bevel blade)
300 µm guarded blade
2 opposed hinges

Ideal for DMEK following Pr Muraine's technique

Available in: 8.00 mm



Shape of the incision

17209D800



KERATOPLASTY Punches

► Deep Well Guarded Punch

Deep Well Guarded Punch
With aspiration
Straight-walled cut (symmetrical double-bevel blade)
Guarded blade

Wide and deep well: ideal for DMEK graft preparation

Available in:


- for peripheral scoring : 9.50 and 10.00 mm
- for diameter sizing : 7.50, 7.75, 8.00, 8.25, 8.50 mm

Full punch : 17212DXXX

Stand-alone blade : 17213DXXX



NEW



Shape of the incision

17212DXXX
17213DXXX


► Guarded Punch

With aspiration
Straight-walled cut (symmetrical double-bevel blade)
300µm guarded blade

Ideal for DMEK graft preparation

Available in:

- for peripheral scoring : 9.50 and 10.00 mm
- for diameter sizing : 7.50, 7.75, 8.00, 8.50 mm



Shape of the incision

17207DXXX



KERATOPLASTY

Trephines

► Adjustable Trephine

Adjustable Trephine (for donor and recipient)
Pre-adjustment of the desired depth, symmetrical double-bevel blade.
Centering guide in contact with the cornea.
Each turn corresponds to a blade descent of 250 µm.
Can be used with artificial chamber #17204.

Available in:


- 6.00 mm
- 6.50 – 9.00 mm (increment of 0.25 mm)
- 9.50, 10.00 mm

Box of 1 unit

PK-DALK



Bague de pré-réglage



Blade profile

17202DXXX

► Non Adjustable Trephine


Simple Trephine (for donor and recipient)
Straight-walled cut: symmetrical double-bevel blade.
Centering guide in contact with the cornea.
Each turn corresponds to a blade descent of 250 µm.
Can be used with artificial chamber #17204.

Available in:

- 7.00 – 9.00 mm (increment of 0.25 mm)
- 9.50 mm

Box of 1 unit

PK



Blade profile

17201DXXX



KERATOPLASTY

Artificial Chambers

► Artificial Chamber

To be used with trephines
#17201DXXX and #17202DXXX
(as a replacement of punch #17200DXXX)

Box of 1 unit




17204

► Artificial Chamber

Can be used with One Use-Plus, One Use
and CBm microkeratomes (#19175, #19155
and #19303).

Box of 1 unit



19182




KERATOPLASTY

Microkeratome Accessories

► Guide-Ring for One Use-Plus

Single-Use Guide-Ring for One Use-Plus motor
(#19175)

Box of 5 units




One Use-Plus
Automated DSAEK

19187

► Guide-Ring for One Use

Single-Use Guide-Ring for One Use turbine (#19155)

Box of 5 units




One Use
DSAEK

19186

► Guide-Ring for CBm

Single-Use Guide-Ring for CBm turbine (#19303)

Box of 5 units



CBm DSAEK-ALTk

19180




KERATOPLASTY

Microkeratome Accessories

► One Use LC Heads for One Use-Plus and One Use

Single-Use Heads for motors One Use-Plus (#19175) and One Use (#19155)
Available in: 110, 130, 150, 200, 250, 300, 350, 400, 450, 500, 550 microns.

Box of 5 units



One Use-Plus
Automated DSAEK


One Use
DSAEK

19184/XXX

► CBSU Heads

Single-Use Heads for CBm turbine (#19303)
Available in: 50, 90, 110, 130, 200, 250, 300, 350 microns.

Box of 5 units



CBm DSAEK

19178/XXX




KERATOPLASTY

Microkeratome Accessories

► Tubing Sets

Single-use Tubing Set for Artificial Chamber

Box of 10 units



19181

Sterile Tubing Set for ACP

Box of 10 units



19192



REFERENCES INDEX

13300	Straight Cannula, for children and adults - 25 G (0.50 mm).....31	17237x10	Capsulorhexis forceps - Curved - 90° pointed tips.....12	17511x10	Utrata Forceps - Straight.....11
13300	Straight Cannula, for children and adults - 25 G (0.50 mm).....45	17242x10	Drysdale Lens Manipulator.....22	17511x5	Utrata Forceps - Straight.....11
17204	Artificial Chamber.....53	17247x10	Plastic Handle.....15	17512x10	Crozafon Forceps - Curved.....11
17210	ROP Kit: Speculum and double-ended Depressor.....38	17249x10	Bubble Level.....37	17512x5	Crozafon Forceps - Curved.....11
19180	Guide-Ring for CBm.....54	17250x1	Mendez Ring.....23, 37	17513x10	Capsulorhexis forceps - Curved - Corydon type.....11
19181	Tubing Sets.....56	17250x10	Mendez Ring.....23, 37	17513x5	Capsulorhexis forceps - Curved - Corydon type.....11
19182	Artificial Chamber.....53	17251x10	Bores Marker.....37	17514x10	Notched Forceps.....14
19186	Guide-Ring for One Use.....54	17252x10	Capsulorhexis forceps - Curved - 45° pointed tips.....12	17514x5	Notched Forceps.....14
19187	Guide-Ring for One Use-Plus.....54	17253x20	Double-ended Depresso.....39	17515x10	Troutman Forceps - Straight.....16, 43
19192	Tubing Sets.....56	17255x1	ROP Kit: Speculum and loop Depressor....39	17515x5	Troutman Forceps - Straight.....16, 43
30601	Micro-sponges.....31	17255x20	ROP Kit: Speculum and loop Depressor....39	17516x10	Paufigue Forceps.....14
30761	Sponcals.....31	17257x10	ROP Kit: Speculum and loop Depressor....38	17516x5	Paufigue Forceps.....14
143BIS	Retrolbulbar Needle - 25 G (0.50 mm).....31	17258x10	Neonatal Wire Speculum.....6, 39	17517x10	Jeweler Forceps - Straight.....17
17200DXXX	Busin Punch.....50	17260X10	Lester IOL Manipulator.....21	17518x10	Corneal Marker.....40, 49
17201DXXX	Non Adjustable Trephine.....52	17261X10	Seibel Lens Manipulator.....21	17518x5	Corneal Marker.....40, 49
17202DXXX	Adjustable Trephine.....52	17300x5	Busin Spatula.....48	17519x10	Corneal Marker.....40, 49
17207DXXX	Guarded Punch.....51	17301x1	23 G Busin Forceps.....48	17519x5	Corneal Marker.....40, 49
17209D800	Muraine Punch.....50	17301x5	23 G Busin Forceps.....48	17520x10	Corneal Marker.....40, 49
17212DXXX	Deep Guarded Punch.....51	17302x5	Price Hook.....23, 48	17520x5	Corneal Marker.....40, 49
17221x10	Bonn Forceps - Straight.....13, 44	17303x5	90° Spatula.....23, 48	17521x10	Corneal Graft Preparation Forceps.....14, 49
17222x10	Tying Forceps - Straight.....15, 42	17500x10	Mc Pherson Tying Forceps - Angled ...16, 43	17521x5	Corneal Graft Preparation Forceps.....14, 49
17223x10	Sinskey IOL Manipulator.....21	17501x10	Tying Forceps - Curved.....16, 43	17522X10	Relax manipulator.....20, 40
17225x10	Tying Forceps - Curved.....15, 42	17502x10	Capsulorhexis forceps - Straight - 90° pointed tips.....10	17522x5	Relax Manipulator.....40
17227x10	Kratz Wire Speculum.....6, 36	17503x10	Capsulorhexis forceps - Curved - 90° pointed tips - Corydon type.....10	17523X10	LASIK manipulator.....20, 41
17228x10	Lester IOL Manipulator.....21	17504x10	Bonn Forceps - Straight.....13, 44	17531x10	Vannas Scissors - Curved.....27
17229x10	Lens Manipulator phaco chopper.....22	17505x10	Bonn Forceps - Angled.....13, 44	17532x10	Castroviejo Scissors - Curved.....27
17230x10	Manipulator.....21	17506x5	Scissors - Sharp.....26, 45	17533x10	Westcott Scissors - Curved.....27
17231x10	Adjustable Speculum.....7	17507x5	Scissors - Blunt.....26, 45	17534x10	Needle Holder - Curved.....17
17232x10	Thornton-Fine Fixation Ring.....22	17508x10	Adjustable Speculum.....7	17534x5	Needle Holder - Curved.....17
17233x10	Mc Pherson Tying Forceps - Angled ...15, 42			17537X10	Wire Speculum.....6, 36



REFERENCES INDEX

17537X5 Wire Speculum6, 36

17539x10 Scleral Marker - 3.5 mm / 4 mm.....36

17539x5 Scleral Marker - 3.5 mm / 4 mm.....36

17541x10 Epithelial Spatula.....41, 49

19178/XXX CBSU Heads55

19184/XXX One Use LC Heads for One Use-*Plus* and
One Use.....55

EasySet34



ALPHABETICAL INDEX

23 G Busin Forceps 17301x1	48
23 G Busin Forceps 17301x5	48
90° Spatula 17303x5	23
90° Spatula 17303x5	48

A

Adjustable Speculum 17231x10	7
Adjustable Speculum 17508x10	7
Adjustable Trephine 17202DXXX	52
Artificial Chamber 17204	53
Artificial Chamber 19182	53

B

Bonn Forceps - Angled 17505x10	13, 44
Bonn Forceps - Straight 17221x10	13, 44
Bonn Forceps - Straight 17504x10	13, 44
Bores Marker 17251x10	37
Bubble Level 17249x10	37
Busin Punch 17200DXXX	50
Busin Spatula 17300x5	48

C

Capsulorhexis forceps - Curved - 90° pointed tips - Corydon type 17503x10	10
Capsulorhexis forceps - Straight - 90° pointed tips 17502x10	10
Capsulorhexis forceps - Curved - Corydon type 17513x10	11
Capsulorhexis forceps - Curved - Corydon type 17513x5	11
Capsulorhexis forceps - Curved - 90° pointed tips 17237x10	12
Capsulorhexis forceps – Curved - 45° pointed tips 17252x10	12
Castroviejo Scissors - Curved 17532x10	27
CBSU Heads 19178/XXX	55
Corneal Graft Preparation Forceps 17521x10	14, 49
Corneal Graft Preparation Forceps 17521x5	14, 49
Corneal Marker 17518x10	40, 49

Corneal Marker 17518x5	40, 49
Corneal Marker 17519x10	40, 49
Corneal Marker 17519x5	40, 49
Corneal Marker 17520x10	40, 49
Corneal Marker 17520x5	40, 49
Crozafor Forceps - Curved 17512x10	11
Crozafor Forceps - Curved 17512x5	11

D

Deep Guarded Punch 17212DXXX	51
Double-ended Depresso 17253x20	39
Drysdale Lens Manipulator 17242x10	22

E

EasySet Toutes les références?	34
Epithelial Spatula 17541x10	41, 49

G

Guarded Punch 17207DXXX	51
Guide-Ring for CBm 19180	54
Guide-Ring for One Use 19186	54
Guide-Ring for One Use-Plus 19187	54

J

Jeweler Forceps - Straight 17517x10	17
--	----

K

Kratz Wire Speculum 17227x10	6, 36
---	-------

L

LASIK manipulator 17523X10	20
LASIK Manipulator 17523x10	41
Lens Manipulator phaco chopper 17229x10	22
Lester IOL Manipulator 17228x10	21
Lester IOL Manipulator 17260X10	21

M

Manipulator 17230x10	21
Mc Pherson Tying Forceps - Angled 17233x10	15, 42
Mc Pherson Tying Forceps - Angled 17500x10	16, 43

Mendez Ring 17250x1	23, 37
Mendez Ring 17250x10	23, 37
Micro-sponges 30601	31
Muraine Punch 17209D800	50

N

Needle Holder - Curved 17534x10	17
Needle Holder - Curved 17534x5	17
Neonatal Wire Speculum 17258x10	6, 39
Non Adjustable Trephine 17201DXXX	52
Notched Forceps 17514x10	14
Notched Forceps 17514x5	14

O

One Use LC Heads for One Use-Plus and One Use 19184/XXX	55
---	----

P

Paufigue Forceps 17516x10	14
Paufigue Forceps 17516x5	14
Plastic Handle 17247x10	15
Price Hook 17302x5	23, 48

R

Relex manipulator 17522X10	20, 40
Relex Manipulator 17522x5	40
Retrobulbar Needle - 25 G (0.50 mm) 143BIS	31
ROP Kit: Speculum and double-ended Depressor 17210	38
ROP Kit: Speculum and loop Depressor 17255x1	39
ROP Kit: Speculum and loop Depressor 17255x20	39
ROP Kit : Speculum and loop Depressor 17257x10	38

S

Scissors - Blunt 17507x5	26, 45
Scissors - Sharp 17506x5	26, 45
Scleral Marker - 3.5 mm / 4 mm 17539x10	36
Scleral Marker - 3.5 mm / 4 mm 17539x5	36
Seibel Lens Manipulator 17261X10	21



ALPHABETICAL INDEX

Sinskey IOL Manipulator [17223x10](#) 21

Sponcals [30761](#) 31

Straight Cannula, for children and adults - 25 G (0.50 mm)
[13300](#)..... 31

Straight Cannula, for children and adults - 25 G (0.50 mm)
[13300](#)..... 45

T

Thornton-Fine Fixation Ring [17232x10](#)..... 22

Troutman Forceps - Straight [17515x10](#)..... 16, 43

Troutman Forceps - Straight [17515x5](#)..... 16, 43

Tubing Sets [19181](#) 56

Tubing Sets [19192](#) 56

Tying Forceps - Curved [17225x10](#)..... 15, 42

Tying Forceps - Curved [17501x10](#)..... 16, 43

Tying Forceps - Straight [17222x10](#)..... 15, 42

U

Utrata Forceps - Straight [17511x10](#) 11

Utrata Forceps - Straight [17511x5](#) 11

V

Vannas Scissors - Curved [17531x10](#) 27

W

Westcott Scissors - Curved [17533x10](#) 27

Wire Speculum [17537X10](#) 6, 36

Wire Speculum [17537X5](#) 6, 36

Fabricant

MORIA SA

15, rue Georges Besse

92160 Antony

FRANCE

Phone: +33 (0) 1 46 74 46 74

Fax: +33 (0) 1 46 74 46 70

www.moria-surgical.com

MORIA Inc

1050 Cross Keys Drive

Doylestown, PA 18902

USA

Phone: (800) 441 1314

Fax: + 1 (215) 230 7670

www.moria-surgical.com

MEDINOVA
> oftalmologia

via Antiniana, 21 | 80078 Pozzuoli | NA (Italy)
tel. +39 081 770 73 27 / +39 081 770 31 98 | fax +39 081 546 32 99
info@medinovasrl.com / www.medinovasrl.com

Moria
Ophthalmic Instruments

Moria Japan K.K.

Arcadia Building 6F

1-12-3 Kanda

SudachoChiyoda-Ku

Tokyo 101-0041

JAPAN

Phone: 81-3-6260-8309

Fax: 81-3-6260-8310

www.moriajapan.com

Moria Commercial (CHINA) CO., LTD.

上海目利亚贸易有限公司

Room H 6F Kaili Building NO.432

West Huai Hai Road, Changning

district, Shanghai, 200052, P.R.C

中国上海市长宁区淮海西路432号

凯利大厦6楼08室

Tel/Fax: +86 021 52586095

www.moria-surgical.com